



# SPERM STAINS



## STAINS FOR THE MORPHOLOGICAL ANALYSIS OF SPERM CELLS

Morphological deviations can substantially affect the fertility of sperm cells. Abnormalities of the head and acrosome in particular impair the fertilization ability of semen cells.

Compared to the so-called “wet mount” method which allows morphological analysis of sperm cells without staining, the use of specific stains provides a better differentiation of the sperm cell regions. Morphological abnormalities can be determined more precisely. For this purpose, two efficient and simple-to-use sperm stains are available from Minitube.

### Spermac Stain

Set for diagnostic staining of sperm cells of all domestic mammals, consisting of

- 50 ml red liquid (Spermac “A”)
- 50 ml pale green liquid (Spermac “B”)
- 50 ml dark green liquid (Spermac “C”)
- 50 ml fixative (clear) liquid



**Spermac Stain**, 4 x 50 ml

REF. : [15405/0000](#)

Spermac stain is used to clearly visualize head, acrosome, equatorial region, centerpiece and tail, so that morphological abnormalities of all regions can be identified.

The head of the sperm cell appears red colored, acrosome, centerpiece and tail are green and the equatorial zone is pale green.

The identification of damage to head and acrosome as well as tail abnormalities are easily detected. Due to the distinctive staining of the acrosome, Spermac is particularly suitable for the identification of the acrosome of equine spermatozoa. Spermac can be used for staining of frozen thawed sperm containing glycerol.

### Instructions for use

1. Prepare a thin smear of diluted semen on a clean slide
2. Fix the smear by dipping it into the fixative for at least 5 minutes
3. Dry on a heating plate at 37 °C for 15 minutes or at ambient temperature over night
4. Staining
  - wash by dipping into tap water
  - dip in stain A for 1 minute and wash by dipping into tap water
  - dip in stain B for 1 minute and wash by dipping into tap water
  - dip in stain C for 1 minute and wash by dipping into tap water
  - dry the smear for approx. 12 hours
5. Microscope analysis with oil immersion (1000 x, phase contrast)




Morphological evaluation of at least 100, preferably 200 sperm cells. Morphologically intact and abnormal sperm cells are counted. The type of abnormality is determined, and depending on their morphological deviation the abnormal cells are divided into different groups:

Acrosome deviations (detached, deformation, ...), head damages (lance form, round, narrow, pear form, ...), neck damages (breaks, plasma droplets), centerpiece and tail deviations (persisting plasma droplets, rolled up, bent, rudimental).

## Accessory

**Dye-bath „Coplin“**

REF. : 15405/2470



Sperm cells, stained with Spermac. The acrosome is clearly visible.

## Farely stain

Staining kit for the diagnostic staining of fresh/diluted (but not frozen) sperm cells of all domestic mammals (including poultry).

Kit consists of 3 components, 250 ml each:

- Bottle A: Fixation Solution
- Bottle B: Aniline-Blue
- Bottle C: Crystal-Violet

**Farely stain, 3 x 250 ml**

REF. : 15405/0026




The blue-violet color contrast allows a differential visualization of the head, the acrosome, the equatorial segment, the centerpiece and the tail of sperm cells. A clear distinction between morphologically intact and abnormal sperm cells is possible as well as the identification of the type of abnormality. Farely is suitable for staining of diluted sperm, unless the extender contains Glycerol. For frozen thawed sperm we recommend Spermac.



### Instructions for use

1. Prepare a thin smear on a clean slide
2. Staining of the air-dried smear
  - dip the smear in stain A for 10 seconds
  - dip in stain B for 20 seconds and wash by dipping into tap water
  - dip in stain C for 5 seconds and wash by dipping into tap water
  - dry the preparation for 12 hours
3. Microscope analysis using oil immersion (1000x, phase contrast)



*Sperm cells of a boar, stained with Farelly*

### Supravital sperm stain with Eosin

The Eosin stain determines the percentage of living and dead cells in semen samples of all domestic mammals. It is simple, rapid and effective. Due to the damaged membrane, dead sperm cells absorb the stain, while living cells remain uncolored. The Eosin stain is not suitable for the evaluation of frozen thawed sperm containing glycerol.

**Eosin G**, 50 ml

REF. : **15405/0025**



### Instructions for use

1. Two drops of Eosin are placed on the edge of a clean, preheated slide
2. One drop of semen is placed beside the Eosin on the slide
3. Carefully mix and prepare the smear
4. Dry on warm tray (approx. 20 minutes)
5. Microscope analysis of 200 cells (at 400x magnification) within the next 30 minutes
6. Calculate the percentage of uncolored, living cells





## Supravital sperm stain with eosin-nigrosin

The combination of Eosin and Nigrosin is used for the supravital stain acc. to Bloom for sperm cells of all domestic mammals.

The staining acc. to Bloom produces a stronger contrast compared to Eosin staining: Dead sperm cells with a damaged plasma membrane are colored by Eosin, while living cells don't absorb the stain. Nigrosin creates a dark background, which simplifies the evaluation. This stain is not suitable for the analysis of frozen thawed sperm containing glycerol.

**Eosin G**, 50 ml REF. : 15405/0025

**Nigrosin**, 50 ml REF. : 15405/0029



### Instructions for use

1. One drop of Eosin, two drops of Nigrosin and one drop of semen are placed on the edge of a clean slide
2. Preparation of smear and microscope analysis similar to the Eosin procedure

