



# PHASE CONTRAST MICROSCOPE FOR SEMEN ANALYSIS



## MBL 2000 Microscope - The universal microscope with binocular or trinocular head, with special configuration for semen analysis

- Excellent phase contrast and bright field optics
- Rugged microscope stand, well engineered mechanics
- All models with built-in heating system and control unit for stage temperature setting
- Available with economic video and monitor system
- Made in Germany

Binocular phase contrast microscope MBL 2000, phase contrast objective 20x, integrated heating system HT 10

REF. : [12010/1000](#)



Binocular phase contrast microscope MBL 2000, phase contrast objectives 10x, 40x, 100x, integrated heating system HT 50

REF. : [12010/2000](#)

Binocular phase contrast microscope MBL 2000, phase contrast objectives 10x, 40x, 100x, integrated heating system HT 50, with ocular color camera and monitor

REF. : [12010/2100](#)

Trinocular phase contrast microscope MBL 2000, phase contrast objectives 10x, 20x, 40x, 100x, integrated heating system HT 200

REF. : [12010/2102](#)

