

TECHNICAL DATA SHEET

REPROJELLY

Product description

ReproJelly is a sperm-safe lubricating gel. It is recommended for Al, rectal palpation and lubrication of AVs in animal species.

Properties

- Composition: purified water, hydroxyethyl cellulose, bacteriostatic and fungicidal agents, stabilizers
- Appearance: clear to light yellow gel
- Odor: neutral
- Produced under GMP norms; components are certified corresponding to Ph EUR
- No ingredients of animal origin
- Iso-osmotic
- Tested for sperm compatibility
- Viscosity optimized to coat gloves, Al equipment and instruments
- Storage: +5 to +25°C

Order information			
ReproJelly, 100 ml tube	REF.	:	11907/0100
ReproJelly, 3 I canister	REF.	:	11907/3000
Dispenser for 3 I canister	REF.	:	11907/3100





