



## TECHNICAL DATA SHEET

# MINITUBE CENTRIFUGE

### Product description

Centrifuge for semen production.  
With swing out rotor and buckets for three different vial sizes.

### Application

Ideal for processing of horse ejaculates.  
Suitable also for blood and faeces work.

### Properties

- Swing out rotor for semen centrifugation
- Maximum RCF: 2200 g
- Speed adjustable between 500 and 4000 rpm (10 rpm steps)
- Memory: 99 programs
- Timer: 0 – 99 minutes & hold (30 sec. steps)
- Weight: 22.8 kg (without rotor)
- Dimensions: 350 x 480 x 340 mm
- Power requirement: 125 W

#### Tube types:

100 ml, round shape, 45 x 100 mm – 4 tubes, Ref. 14602/0801

50 ml, point shape, 29 x 115 mm – 4 tubes, Ref. 14607/0053

15 ml, point shape, 17 x 120 mm – 8 tubes, Ref. 17500/0201

### Order information

**Centrifuge** REF. : 14601/0000



minitube